



L7 INFORMATICS



Release Notes for L7|ESP®

Version 3.4.0.0

Author: [Matthew Garwin](#)

November 2024

Table of contents

New Features.....	3
Custom Fields.....	3
Vector Fields in L7 MES.....	3
L7 Notebooks.....	3
Naming Notebook Widgets.....	3
Autosave Entries to avoid data loss.....	4
Audit History in Notebooks.....	5
Increased Visualizations.....	6
PDF generation for Notebook Entries.....	6
Rich Text Editor.....	7
Signatures in Notebook Entries.....	8
Extending widget features for Notebook entry.....	10
L7 HUB.....	12
Support (Bundle) approvals across systems.....	12
Support for Dynamic Execution Plans.....	13
L7 MASTER.....	14
General Dynamic Execution Plans.....	14
BOM (Bill of Materials) Dynamic Execution Plan.....	15
Dynamic Execution Plans - Save As.....	15
Vector Fields.....	15
Protocol Actions.....	15
L7 MES.....	16
Looping and Repeating Batch Sections.....	16
Vector Field Support.....	18
Improved Tablet Support in EBR.....	18
Applying Dynamic Execution Plans in MES.....	18
L7 Scheduling.....	19
Archive Availability.....	19
Archive Tasks in the Calendar Tab.....	19
Create and Maintain Auto-generated Tasks.....	20
Maintain and Update Resource Availability.....	21
Resource Management Tool in Monitor.....	21
Experiments in L7 LIMS.....	22
Applying Dynamic Execution Plans when creating an Experiment for L7 LIMS.....	22
Resolved Issues.....	23
Known Issues.....	23
Constraints.....	24
Seeding Users.....	24

New Features

Custom Fields


Vector Fields in L7 MES

Introduced functionality to allow for tables of data, with columns that can be set as text or numbers. Vector fields can also be configured as time series, with a date/timestamp as the first column, if desired.

L7 Notebooks

Naming Notebook Widgets

L7 Notebooks now provides the ability to auto-generate unique names for Advanced Table, Widget Container and Spreadsheet widgets when these are added to the entry and/or template. Additionally, these auto-generated names can be modified as needed to provide more meaningful names to the widgets. To make it more flexible and intuitive, the application warns the user if there are duplicate or empty names for the widgets but does not block the user from saving such entries.



The screenshot displays the L7 Notebooks interface. At the top, there is a navigation bar with 'NOTEBOOKS', 'My Entries', 'Library', 'Entries', 'Templates', and 'Search'. Below this, the current entry is identified as 'Lib-1 > NB Entry 001'. The interface is divided into two main sections: 'Entry Details' on the left and the 'Editor' on the right.

The 'Entry Details' sidebar includes fields for Name (NB Entry 001), Description, Tags (with a placeholder 'Enter tags here'), Signature Flow (None), Owner (system admin), Created (1 minute ago), Updated (1 minute ago), Access (Owner), and Entry Status (Edited by system admin). The Editor section has tabs for 'Editor', 'Report', and 'History'. It contains two tables:

	A	B	C
1			

	A	B	C	D
1				

Below the tables are expandable sections: 'Widget container 1', 'Widget container 2', and 'Sheet1'. The 'Sheet1' section is expanded to show a spreadsheet widget with a menu (Home, Insert, Formulas, View, Export) and a grid of cells (A1 to P1).

An orange warning icon in the bottom right corner of the spreadsheet widget indicates an error: 'Spreadsheet Widget: Spreadsheet widget names should not be empty'.

At the bottom of the interface, there is a footer with the version '3.4.0.0.20241018114208', a development notice 'DEVELOPMENT BUILD. NOT FOR PRODUCTION USE', and a copyright notice 'Copyright © 2013-2024 L7 Informatics. All Rights Reserved'.

Autosave Entries to avoid data loss

L7 Notebooks provides the option to auto-save entries at user defined intervals so that data is not lost in case of system timeout or user navigating out from the entry unintentionally. The auto-save can be enabled, and the interval can be configured through the config app.

Lib-1 > NB Entry 001

Review | Archive | Request for Review | Save

Last Auto Saved: a few seconds ago

Entry Details: Name: NB Entry 001, Description: [Empty], Tags: [Enter tags here], Signature Flow: None, Owner: system admin, Created: 20 minutes ago, Updated: 1 minute ago, Access: Owner, Entry Status: Edited by system admin

Editor | Report | History

Show Header

Entry Auto-saved and auto save message displayed at the top right of this entry page. Autosave can be set on/off from the ESP configuration.

Table 1			
	A	B	C
1			

Table-New				
	A	B	C	D
1				

Widget container 1

There are no widgets yet

Drag and drop widget from right panel here to add it to the content

Widgets Palette: Header, Subheader, Paragraph, Bullet list, Ordered list, Checklist list, Table, Spreadsheet, Rich text, Generic renderer, Horizontal separator, Date separator, Image, File, Widget container

3.4.0.0.20241018114208 DEVELOPMENT BUILD. NOT FOR PRODUCTION USE Copyright © 2013-2024 L7 Informatics. All Rights Reserved

Audit History in Notebooks

L7 Notebooks now provides the capability to select a printer and physically print the audit trail (records in the History page) for the Library, Entries, Entry Details, Templates and Template Details pages.

Entry Details: Name: NB Entry 001, Description: [Empty], Tags: [Enter tags here], Signature Flow: None, Owner: system admin, Created: 20 minutes ago, Updated: 1 minute ago, Access: Owner, Entry Status: Edited by system admin

Print | 1 page | Last Auto Saved: 3 minutes ago

Destination: Save as PDF

Pages: All

Layout: Portrait

More settings

Cancel | Save

Details

Details

Details

11/20/24, 4:29 PM NB Entry 001 history

Date, Action, Performed By

2024-10-30T04:03:48.418078Z, created Notebook Entry "NB Entry 001" (bde92cf8-1289-42cd-8235-8fc2aac3b361), admin@localhost

2024-10-30T04:06:28.587237Z, created Spreadsheet widget "Sheet1", admin@localhost

2024-10-30T04:06:28.587237Z, created Spreadsheet widget "Table 1", admin@localhost

2024-10-30T04:06:28.587237Z, created Table widget "Table 2", admin@localhost

2024-10-30T04:06:28.587237Z, created Table widget "Table 1", admin@localhost

2024-10-30T04:06:28.587237Z, created Spreadsheet widget "Sheet2", admin@localhost

2024-10-30T04:20:48.437239Z, updated Spreadsheet widget name from "Sheet2" to "(empty)", admin@localhost

2024-10-30T04:20:48.437239Z, updated Table widget name from "Table 2" to "Table-New", admin@localhost

2024-10-30T04:24:28.378669Z, created Header widget "Entry Auto-saved and...", admin@localhost

2024-10-30T04:25:47.285163Z, updated Header widget "Entry Auto-saved and..." content, admin@localhost

3.4.0.0.20241018114208 DEVELOPMENT BUILD. NOT FOR PRODUCTION USE Copyright © 2013-2024 L7 Informatics. All Rights Reserved

Increased Visualizations

L7 Notebooks now supports the feature to create charts using data from the spreadsheet widget.

The screenshot displays the L7 Notebooks interface. On the left, the 'Entry Details' sidebar shows information for 'NB Entry 002', including its name, description, tags, signature flow, owner (system admin), creation and update times (1 minute ago), access level (Owner), and entry status (Edited by system admin). The main area features a 'Spreadsheet' widget with a table containing the following data:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	1	100	1000													
2	2	200	2000													
3	3	300	3000													
4	4	400	4000													
5	5	500	5000													
6																
7																
8																
9																
10																
11																
12																
13																

Below the spreadsheet is a 'Graph-1' widget, which is a line chart titled 'Sheet1 - A x B'. The chart displays a linear trend with the following data points:

A	B
1	100
2	200
3	300
4	400
5	500

The interface also includes a top navigation bar with 'NOTEBOOKS', 'My Entries', 'Library', 'Entries', 'Templates', and 'Search'. The bottom status bar shows the version '3.4.0.0.20241018114208', a development build notice, and the copyright 'Copyright © 2013-2024 L7 Informatics. All Rights Reserved'.

PDF generation for Notebook Entries

Expanded capability in L7 Notebooks to download a .zip file which contains a PDF version of the latest notebook entry version including all attached images and files in their native formats.

The screenshot displays the L7 Notebooks interface for an entry titled "NB Entry 002". The interface includes a top navigation bar with "NOTEBOOKS", "My Entries", "Library", "Entries", "Templates", and "Search". The main content area is divided into three sections: "Entry Details", "Report", and "History".

Entry Details:

- Name: NB Entry 002
- Description: (Empty text area)
- Tags: Enter tags here
- Signature Flow: None
- Owner: system admin
- Created: 12 minutes ago
- Updated: 1 minute ago
- Access: Owner
- Entry Status: Edited by system admin

Report Section:

PDF generated 10/30/2024 04:42 am UTC

Open by system admin 10/30/2024 04:29 am UTC

Table:

Sheet1			
1	100	1000	
2	200	2000	
3			
4			
5			

Graph-1:

A line graph showing a linear trend. The x-axis is labeled 'A' and ranges from 1 to 5. The y-axis ranges from 0 to 600. The data points are:

A	Value
1	100
2	200
3	300
4	400
5	500

Exporting Entry Dialog:

Exporting Entry will generate .zip file with latest version of Entry, the latest version of attached file(s) and download to your local computer. This may take a while. While exporting is in progress, you may continue to use other areas of ESP. Would you like to continue?

Buttons: Cancel, Continue Export


Export Settings:

- Entry:
 - Paper Size: Letter
 - Pages: All
 - Margin: Default (20mm)
- Worksheet/Table:
 - Name: Select Worksheet or Table
 - Scale: (Dropdown menu)

Footer: 3.4.0.0.20241018114208 DEVELOPMENT BUILD. NOT FOR PRODUCTION USE. Copyright © 2013-2024 L7 Informatics. All Rights Reserved.

Rich Text Editor

L7 Notebooks now supports the capability to add nested bullets and number lists in its Rich Text Editor widget. Additional configuration has been added to hide or show the Table Creation option within the rich text editor (Table Creation icon is set to be hidden by default).



Signatures in Notebook Entries

Users can now sign Notebook Entries asynchronously. Email notifications for asynchronous signatures will all be sent at the start of the signature flow while Email notifications for synchronous signatures will continue to be sent sequentially.

L7 Notebooks has added the ability to sign in as another signer in cases when a specific signer is unavailable to sign the entry or is acting as a witness to approving details in the experiment. This feature is driven through the configuration flag, `other_user_approval` which is disabled by default. If the flag is enabled, a user can sign an entry as another user who is part of the currently assigned signature flow.

Notebook Entry Review: NB Entry 002
Please review the following information

Select an action: **Approve** | Reject

Sign as: 1. Owner: system admin

Reason:

Username: 1. system admin (Owner) | 2. user1 (User)

Password:

Comment:

I understand that my approval will be an electronic representation of my signature for all purposes, similar to that of a pen-and-paper signature or initial.

Buttons: Cancel | Approve

Background Entry Details:
Name: NB Entry 002
Description:
Tags:
Owner: system admin
Created: 10/29/2024 09:29 pm PDT
Last Updated: 10/29/2024 09:53 pm PDT
Entry Status: Under Review by system admin

Background Spreadsheet (Sheet1):

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R
1	1	100	1000															
2	2	200	2000															
3	3	300	3000															
4	4	400	4000															
5	5	500	5000															

To avoid unnecessary edits to the entry during sign-off, a capability has been added to block signers from unlocking the entry under review. The capability to unlock a notebook entry under review is now limited to only the Owner and Admins of the entry.

Background Entry Details:
Name: NB Entry 002
Description:
Tags:
Owner: system admin
Created: 25 minutes ago
Updated: 1 minute ago
Access: Can Edit
Entry Status: Under Review by system admin


Background Spreadsheet (Sheet1):

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R
1	1	100	1000															
2	2	200	2000															
3	3	300	3000															
4	4	400	4000															
5	5	500	5000															

Extending widget features for Notebook entry

Added functionality to expand/collapse widgets for better editing and better readability of Notebook entries.

Added the capability to hide or show a preview of an uploaded file to improve the editing experience.



The screenshot displays a web-based interface for editing a notebook entry. The top navigation bar includes 'NOTEBOOKS', 'My Entries', 'Library', 'Entries', 'Templates', and 'Search'. The user is logged in as 'esp | system admin'. The entry is titled 'Lib-1 > NB Entry 001'. The left sidebar shows 'Entry Details' with fields for Name, Description, Tags, Signature Flow, Owner, Created, Updated, Access, and Entry Status. The main editor area has tabs for 'Editor', 'Report', and 'History'. It shows two widget containers, each with a collapse/expand arrow. The first container is expanded to show a preview of a CSV file named 'Sheet90.csv'. The CSV data is as follows:

#	A	B	C	D	E
1	header-1	header-2	header-3	header-4	header-5
2	1	10	100	1000	10000
3	2	20	200	2000	20000
4	3	30	300	3000	30000
5	4	40	400	4000	40000
6	5	50	500	5000	50000
7	6	60	600	6000	60000
8	7	70	700	7000	70000
9	8	80	800	8000	80000
10	9	90	900	9000	90000
11	10	100	1000	10000	100000
12					
13					
14					
15					
16					
17					

Below the CSV preview, there is a collapsed widget for 'Book4.xlsx' and another for 'Sheet1'. The footer of the interface shows the version '3.4.0.0.20241018114208', a development build notice, and the copyright '© 2013-2024 L7 Informatics, All Rights Reserved'.

Users can now create nested ordered lists and bullet lists within the notebook entry.

The screenshot shows the L7 Informatics interface for a notebook entry titled "NB Entry 001". The interface includes a top navigation bar with "NOTEBOOKS", "My Entries", "Library", "Entries", "Templates", and "Search". The main content area is divided into three sections: "Entry Details", "Editor", and "Widgets Palette".

Entry Details: Name: NB Entry 001, Description: (empty), Tags: (empty), Signature Flow: None, Owner: system admin, Created: 1 minute ago, Updated: 1 minute ago, Access: Owner, Entry Status: Edited by system admin.

Editor: Shows a list of items and a numbered list. The list of items includes Item1, Item2, Item3, Item4, Item5, Item6, and Item7. The numbered list includes Test-1 through Test-10.

Widgets Palette: Lists various widgets such as Header, Subheader, Paragraph, Bullet list, Ordered list, Checklist, Table, Spreadsheet, Rich text, Generic renderer, Horizontal separator, Date separator, Image, File, and Widget container.

Introduced functionality to duplicate the advanced table widget along with the format and data in the table widget and use the duplicated widget as a new independent table widget.

The screenshot shows the L7 Informatics interface for a notebook entry titled "NB Entry 002". The interface includes a top navigation bar with "NOTEBOOKS", "My Entries", "Library", "Entries", "Templates", and "Search". The main content area is divided into three sections: "Entry Details", "Editor", and "Widgets Palette".

Entry Details: Name: NB Entry 002, Description: (empty), Tags: (empty), Signature Flow: None, Owner: system admin, Created: 1 minute ago, Updated: 1 minute ago, Access: Owner, Entry Status: Edited by system admin.


Editor: Shows a table widget and its duplicate. The table has columns A, B, and C, and rows 1 and 2. The data is as follows:

	A	B	C
1	New-1	Item-1	
2	New-2	Item-2	

The duplicate is titled "Copy of Table 1" and contains the same data as the original table.

Widgets Palette: Lists various widgets such as Header, Subheader, Paragraph, Bullet list, Ordered list, Checklist, Table, Spreadsheet, Rich text, Generic renderer, Horizontal separator, Date separator, Image, File, and Widget container.

Added ability to expand the spreadsheet widget offering a better view to edit data in the widget.



L7|HUB®

Support (Bundle) approvals across systems

Administrators now have the ability to see and filter on the review status of remote bundles. Installation of remote bundles may be based on review status, based on configuration.

Support for Dynamic Execution Plans

Dynamic Execution Plans may now be included in collections and bundles, and can be transferred between systems through L7|HUB.

General Dynamic Execution Plans

L7|ESP now supports General Execution Plans, which enable run-time modifications to the configuration of any fields in a Workflow Chain.

The screenshot shows the 'Dynamic Execution Plans: LML - General US' interface. The main area displays a tree view of the workflow chain with the following structure:

- Thaw Large Molecule
 - Thaw

Display Name	Type	Default Value	Required
Donor Tracking Number - US	Free-form text entry	-	Yes
Number Vials - General	Numbers only allowed	-	Yes
Thaw Time - General	Numbers only allowed	-	No
 - 5L Bioreactor Preparation
 - 5L Bioreactor Prep

Display Name	Type	Default Value	Required
Number Vials Added	Numbers only allowed	-	Yes
Initial Jacket Temperature	Numbers only allowed	-	No
 - 5L Bioreactor Monitoring
 - Bioreactor Monitoring

Display Name	Type	Default Value	Required
Readings Instructions	Instructions	-	No
Operator Readings	Vector	-	No
 - 10L Bioreactor Transfer
 - 10L Bioreactor Transfer

Display Name	Type	Default Value	Required
Date/Time of Transfer	Date/Time	-	Yes
 - 10L Bioreactor Monitoring
 - Bioreactor Monitoring

Display Name	Type	Default Value	Required
Viable Cell Density	Sample Point	-	No

The right-hand side of the interface shows the 'Inspector: Readings Instructions' panel with various configuration options like 'Basic Details', 'ID', 'Tags', 'Type', 'Grouping', 'Instructions', 'Shared', 'Visible', 'Read-only', 'Required', 'Reportable', 'Parameters', and 'Event Handler'. The 'Instructions' field contains the following text:

```
<p>Record the Jacket Pressure and Temperature, measure the pH every three hours.</p><p style="color: rgb(255, 0, 0);">If the pH is too low, add .11 L of Buffer.</p>
```


This screenshot shows the same interface as above, but with a text editor window open over the 'Readings Instructions' field. The text editor contains the following text:

```
Record the Jacket Pressure and Temperature, measure the pH every three hours.
If the pH is too low, add .11 L of Buffer.
```

The text editor includes a rich text toolbar with options for bold, italic, underline, strikethrough, text color, background color, bulleted list, numbered list, link, unlink, indent, outdent, and other standard text formatting tools. The status bar at the bottom of the editor shows 'CHARS: 98 WORDS: 23 POWERED BY JOOT'.

BOM (Bill of Materials) Dynamic Execution Plan

L7|ESP now supports BOM (Bill of Materials) Execution Plans, which enable run-time modifications to the configuration of Inventory Item Use fields and the Resource Link to Container fields that are referenced by them.



Dynamic Execution Plans - Save As

Users may now save a copy of a Dynamic Execution Plan in order to create a similar Execution Plan more easily.

Vector Fields

Vector fields can be added to protocols intended for use in L7 MES (Use in Manufacturing = Yes). These fields are structured as tables, with columns that can be set as text or numbers. Vector fields can also be configured as time series, with a date/timestamp as the first column, if desired.

Numeric columns can be configured with a specification type, units of measure, decimal digits, and minimum/maximum ranges.

The vector table can be configured with a minimum, maximum, and default number of rows.


Protocol Actions

The Protocol Actions **Create MES Batch** and **Start MES Workflow Chain** now support all four Dynamic Execution Plans types: General, BOM, Specification, and Sample.

L7 MES

Looping and Repeating Batch Sections

L7 MES Recipes (Workflow Chains) can be configured to re-execute Batch Sections (Workflows), by repeating the same Batch Section one or more times, or by looping backward in the Chain. When a Batch Section is re-executed, the system appends a (sequential) letter to the Section's number.



Step 3B.1: Bioreactor Monitoring

Step 3B.1.1: Readings

Record the Jacket Pressure and Temperature, measure the pH every two hours.
If the pH is too low, add 1 L of Buffer.

Operator Readings

Date/Time Recorded	Jacket Temperature (°C) [78...84]	Pressure [20...23]	pH [6.8...7.6]	Amt Buffer Added (L) [05...15]	Operator
2024-11-02 06:15 CDT	79.9	21.4	▲ 7.7	0.08	boperator
2024-11-02 08:15 CDT	80.0	21.6	7.3	N/A	boperator
2024-11-02 10:14 CDT	80.0	21.6	7.4	N/A	boperator
2024-11-02 12:15 CDT	80.0	21.7	7.6	N/A	boperator
2024-11-02 14:16 CDT	80.1	21.9	7.5	N/A	boperator
Select date and time					

At 20 hours, draw 3 samples and send to the lab for Viable Cell Density.

Viable Cell Density Samples: 3

Samples Taken	Time Registered
ESP06724	10/31/2024 10:59 am CDT
ESP06725	10/31/2024 10:59 am CDT
ESP06726	10/31/2024 10:59 am CDT

Users may sort the data displayed in all tabs (except EBR) by numeric order or by order of execution. All reports also observe the tab-based sort order.

MANUFACTURING EXECUTION SYSTEM

Batch Number: LMS-0031 | EBR | Sample Plan | BOM | Specifications | Quality Review | Sign Off | Reports | Attachments | History

Specification Plan: LML - Specs US
Batch Section: All

Production Record: Large Molecule - Looping

Material Number: ALM Capybara
Version: 12.4 (10/31/2024 10:37 am CDT)
Reagent: Tank 8

Field	Spec Type	UOM	Min	Max	Actual Value	Comment
3.1.1 - Operator Readings						
Pressure						
2024-10-31 06:15 UTC	KPP	-	20	23	22.0	-
2024-10-31 08:05 UTC	KPP	-	20	23	24.0	▲
2024-10-31 10:16 UTC	KPP	-	20	23	23.5	▲
2024-10-31 12:46 UTC	KPP	-	20	23	22.9	-
2024-10-31 14:47 UTC	KPP	-	20	23	22.5	-
pH						
2024-10-31 06:15 UTC	KPP	-	6.8	7.6	6.9	-
2024-10-31 08:05 UTC	KPP	-	6.8	7.6	7.5	-
2024-10-31 10:16 UTC	KPP	-	6.8	7.6	7.8	▲
2024-10-31 12:46 UTC	KPP	-	6.8	7.6	7.4	-
2024-10-31 14:47 UTC	KPP	-	6.8	7.6	7.1	-
Amt Buffer Added						
2024-10-31 06:15 UTC	KPP	L	.05	.15	N/A	-
2024-10-31 08:05 UTC	KPP	L	.05	.15	N/A	-
2024-10-31 10:16 UTC	KPP	L	.05	.15	0.10	-
2024-10-31 12:46 UTC	KPP	L	.05	.15	N/A	-
2024-10-31 14:47 UTC	KPP	L	.05	.15	N/A	-
3A.1.1 - Operator Readings						
Jacket Temperature						
2024-11-01 18:06 UTC	CPP	°C	78	84	84.1	▲
2024-11-01 08:12 UTC	CPP	°C	78	84	84.0	-
2024-11-01 10:53 UTC	CPP	°C	78	84	83.9	-
Pressure						
2024-11-01 18:06 UTC	KPP	-	20	23	22.9	-
2024-11-01 08:12 UTC	KPP	-	20	23	23.0	-
2024-11-01 10:53 UTC	KPP	-	20	23	22.6	-

non-versioned release | DEVELOPMENT BUILD: NOT FOR PRODUCTION USE | Copyright © 2019-2024 L7 Informatics. All Rights Reserved.

A new Transition Genealogy Report lists the order of transition between Batch Sections.

Material Number: ALM Capybara
Process: Large Molecule - Looping
Version: 12.4 (10/31/2024 10:37 am CDT)
Batch Number: LMS-0031

L7 INFORMATICS

Transition Genealogy Report

Execution Order

Section	Next Section
1.0 Thaw Large Molecule	2.0 5L Bioreactor Preparation
2.0 5L Bioreactor Preparation	3.0 5L Bioreactor Monitoring
3.0 5L Bioreactor Monitoring	3A.0 5L Bioreactor Monitoring
3A.0 5L Bioreactor Monitoring	3B.0 5L Bioreactor Monitoring
3B.0 5L Bioreactor Monitoring	4.0 10L Bioreactor Transfer
4.0 10L Bioreactor Transfer	5.0 10L Bioreactor Monitoring
5.0 10L Bioreactor Monitoring	5A.0 10L Bioreactor Monitoring
5A.0 10L Bioreactor Monitoring	5B.0 10L Bioreactor Monitoring
5B.0 10L Bioreactor Monitoring	5C.0 10L Bioreactor Monitoring
5C.0 10L Bioreactor Monitoring	5D.0 10L Bioreactor Monitoring
5D.0 10L Bioreactor Monitoring	6.0 Purification 1
6.0 Purification 1	7.0 Purification 2
7.0 Purification 2	Process Finished

Powered by L7 ESP Enterprise Science Platform

Transition Genealogy Report
Printed By: aoper
Printed Date: 10/31/2024 01:35 pm CDT

General Plan: LML - General US
Specification Plan: LML - Specs US
Sample Plan: Default
BOM Plan: Default

Page 1 of 2

Vector Field Support

L7 MES now supports Vector field types.

MANUFACTURING EXECUTION SYSTEM

Section 5.0 10L Bioreactor Monitoring | EBR | Sample Plan | BOM | Specifications | Quality Review | Sign Off | Reports | Attachments | History

General Plan: LML - General US

Production Record: Large Molecule - Looping

Batch Number: LMS-0031 | Material Number: ALM Copybara
Batch Status: in process | Version: 12.4 (10/31/2024 10:37 am CDT)
Patient: PT0004 | Reagent: Tank 8

Step 5.1: Bioreactor Monitoring

Step 5.1.1: Readings

Record the Jacket Pressure and Temperature, measure the pH every two hours.
If the pH is too low, add .1 L of Buffer.

Operator Readings

Date/Time Recorded	Jacket Temperature (°C) (78...84)	Pressure (20...23)	pH (6.8...7.6)	Amt Buffer Added (L) (.05...15)	Operator
2024-11-02 18:01 CDT	79.0	22.0	▲ 7.8	0.10	boperator
2024-11-02 20:05 CDT	79.5	22.4	7.6	N/A	boperator
2024-11-02 22:02 CDT	80.1	22.7	7.4	N/A	boperator
2024-11-03 00:02 CDT	▲ 85.0	22.9	7.4	N/A	asper
2024-11-03 02:03 CST	84.0	21.0	▲ 7.9	0.11	asper
2024-11-03 04:04 CST	83.8	21.4			

Comments

Jacket Temperature, Row 4

Alan Operator
This bioreactor is nearing its recalibration date and may need to be taken out of commission until recalibrated

At 20 hours, draw 3 samples and send to the lab for Viable Cell Density.

Viable Cell Density | Samples: 0

Sample Type: Viable Cell Density (VCD) [3..3] (0)

In EBR, Quality Review, and the final Batch Report, vector fields appear as tables. Users can add and remove rows within the minimum/maximum number of rows defined for the vector table.

Users may add comments to individual cells in the vector table, and may set individual cells as N/A.

Numeric columns are recognized as specifications. Out of range errors are highlighted in EBR and Quality Review, and are included in the Specifications tab and Specifications report.

In Quality Review and the Exceptions Report, exceptions in individual table cells are highlighted.

Improved Tablet Support in EBR

EBR has been modified for an improved tablet experience. Fixed modal windows are now scrollable to allow easier access to buttons while the virtual keyboard is displayed, and other issues with excessive scrolling in EBR have been addressed.

Applying Dynamic Execution Plans in MES


L7 MES now supports four types of Execution Plan: General, BOM, Specification, and Sample. All four may be applied when creating an MES batch. For more information about the new General and BOM Plan types, see the L7|MASTER section of these release notes.

It is possible for the same field to be included in both the General and one of the three specialized Plans. If both are applied, the system will apply all configuration settings that are unique within the two Plan types; where there is overlap, the system will apply the configuration setting from the specialized Plan type.

L7 Scheduling


Archive Availability

When archiving a scheduled resource, users will now be asked if they would like the resource removed from all future scheduled assignments as well.




Archive Tasks in the Calendar Tab

Personnel will be notified via an orange toast notification and an email notification when a task they were assigned to has been archived and the option to send the notification is selected during task archival.




Create and Maintain Auto-generated Tasks

Improvements were made to the Save Confirmation modals along with the notifications for editing Planned Action Chains through the Task Edit modal.



Maintain and Update Resource Availability

L7 Scheduling now supports editing the start date for a future schedule.



Resource Management Tool in Monitor

Streamlined the Inspector Panel in the Resources tab to only expand the details for the specific item selected.

The screenshot shows the 'Monitor' interface for 'Monitor Task Progress'. The main view is a 'Timeline' for 'Experiment-2' (Oct 29 - 30). A task 'WF1' is shown on Tuesday, October 29, 2024. A tooltip for this task provides the following details:

- Task Name: WF1
- Start date: 10/29/2024 10:48 pm
- End date: 10/30/2024 12:03 am
- Duration: 1 hour 15 minutes
- Due date: 10/31/2024

The right-hand panel, titled 'Inspector: Experiment-2 - WF1', displays the following details:

- Experiment: Experiment-2 - WF1
- Description: None
- Tags: None
- Workflow Name: WF1
- Start Time: 10/29/2024 10:48 pm
- End Time: 10/30/2024 12:03 am
- Due date: 10/31/2024
- Personnel: Unassigned
- Workgroup: Unassigned
- Location: Unassigned
- Equipment: Unassigned
- Baseline Start: 10/31/2024 03:30 pm
- Baseline End: 10/31/2024 04:45 pm
- Protocol: P1

Experiments in L7 LIMS

Applying Dynamic Execution Plans when creating an Experiment for L7 LIMS

L7 LIMS now supports four types of Execution Plan: General, BOM, Specification, and Sample. All four may be applied when creating an experiment for L7 LIMS. For more information about the new General and BOM Plan types, see the L7|MASTER section of these release notes.

It is possible for the same field to be included in both the General and one of the three specialized Plans. If both are applied, the system will apply all configuration settings that are unique within the two Plan types; where there is overlap, the system will apply the configuration setting from the specialized Plan type.

Resolved Issues

Known Issues from previous releases resolved in 3.4.0.0

Key	Summary	Application	ZenDesk ID
APPS-19206	Entities App: Custom views not saving/appearing correctly when custom view is set to default view	Entities	N/A
APPS-18986	LIMS: Unable to scan sample barcodes into worksheets (Context lost)	LIMS	N/A

Known Issues

Key	Summary	Application
APPS-19632	DEP creation - WFC should not be a live link	L7 MASTER
APPS-19628	3.3.0.x hotfix CLONE - LIMS left sidebar Filtering throws error	L7 LIMS
APPS-19625	DEP - Save As sometimes saves additional information in original DEP	L7 MASTER
APPS-19618	L7 MASTER > DEP: In the DEP details page, unable to add Workgroups to the DEP (only Default option 'Any' is displayed)	L7 MASTER
APPS-19617	L7 MASTER > DEP: Workgroup pick list is displayed within a very small scrollable section making it difficult to select workgroup options in the DEP creation modal	L7 MASTER
APPS-19599	L7 LIMS left sidebar Filtering throws error	L7 LIMS
APPS-19561	L7 MASTER: Container Type - Unable to remove resource types from Container Type field 'Contains'	L7 MASTER
APPS-19559	L7 MES - Exception Report doesn't include vector field cells with out of range or incomplete	L7 MES
APPS-19556	L7 MES - column headings cut off on some tabs	L7 MES
APPS-19527	L7 MASTER - DEP - Sample point type is displaying in BOM and Spec Plan DEP	L7 MASTER
APPS-19526	L7 MES: DEP not applied to Protocol in 2nd WF	L7 MES
APPS-19517	Various onRender and onChange issues with grouped protocols	L7 LIMS
APPS-19508	L7 MES/Quality Review - Next Step should be read only in Editor Portal	L7 MES
APPS-19449	L7 MES- Vector Columns name updates to match Figma	L7 MES
APPS-19441	L7 MASTER > DEPs: Notify groups unable to be modified	L7 MASTER
APPS-19415	L7 Notebooks: The Auto-save keeps triggering even if the user is not actively working on any entry	L7 Notebooks
APPS-19351	L7 LIMS: OnRender called twice when switching from Rows/Transpose to Flexview view	L7 LIMS

APPS-19097	L7 MES: No error toast upon creation of EBR with DEP when process WFC had field removed after creation of that DEP which included that field	L7 MES
APPS-18951	L7 MES: EBR > BOM tab > Available Items & Item Usage Tables - Lot expiration dates 2 days before Inventory app expiration dates	L7 MES
APPS-18930	L7 MES: Custom field in Entity Type with DEP Resource Link restricted to a plan & WFC - Field name is cut off in new EBR modal	L7 MES
APPS-16403	Uploading file names with '#' breaks viewing/download ability for that file.	Data
APPS-15918	Duplicate toast notifications when POST is used	L7 ESP Core

Constraints

In order to meet regulatory requirements for end-to-end encryption, L7 Informatics strongly recommends customers run L7|ESP with SSL enabled even in network environments where L7|ESP is not accessible from the internet. This ensures data integrity from the point of data origin through data storage.

L7|ESP is validated only in SSL-enabled environments, running with SSL disabled constitutes running in an unvalidated state and L7 Informatics cannot guarantee that L7|ESP will function properly without SSL enabled in future releases.

Seeding Users

Seeding user files will recognize the workgroups. When seeding a customer repo, workgroups in the user file should be removed. L7|ESP client re-organizes the resources, where users are seeded before the workgroups.

L7|ESP Testing activities for L7|ESP 3.4.0.0 were successfully concluded at 17-Oct-2024.